

BRK™ Transseptal Needles


Product Highlights

- Available in adult sizes
- Bevel angle 50 degrees
- Pointer on the needle shield shows the direction of the needle curve
- Specifically designed to work with Abbott fixed curve and steerable transseptal introducers (except model G408332)

Ordering Information

Stainless Steel Transseptal Needle and Stylet (1 unit per box)

REORDER NUMBER	CURVE TYPE	SIZE	NEEDLE USABLE LENGTH (CM)
407200	BRK™	Adult, 18 ga.	71
407201	BRK-1™	Adult, 18 ga.	71
407205	BRK™	Adult, 18 ga.	89
G407215	BRK-1™	Adult, 18 ga.	89
407206	BRK™	Adult, 18 ga.	98
407207	BRK-1™	Adult, 18 ga.	98



BRK™ Transseptal Needles, XS Series

Product Highlights

- Bevel angle 30 degrees
- Steeper primary bevel angle and two back bevels that combine to form a distinct point at the tip of the needle
- Pointer on the needle shield shows the direction of the needle curve
- Specifically designed to work with Abbott fixed curve and steerable transseptal introducers (except model G408332)

Ordering Information

Stainless Steel Transseptal Needle and Stylet (1 unit per box)

REORDER NUMBER	CURVE TYPE	SIZE	NEEDLE USABLE LENGTH (CM)
G407208	BRK™ XS	Adult, 18 ga	71
G407209	BRK-1™ XS	Adult, 18 ga	71
G407210	BRK™ XS	Adult, 18 ga	89
G407216	BRK-1™ XS	Adult, 18 ga	89
G407211	BRK™ XS	Adult, 18 ga	98
G407212	BRK-1™ XS	Adult, 18 ga	98

