

# GUIDERIGHT™

0.035" Guide Wire for Accurate Access

## GUIDING THROUGH ACCESS

GuideRight™ is a stainless steel access Guide Wire, designed to overcome standard clinical challenges

### **Stainless steel core**

- Ideal for general needle access, introduction of sheaths, and catheter placement

### **PTFE (Teflon) pre-coating**

- Exceptional surface properties with increased slide capability and reduced friction between the catheter and Guide Wire

### **Versatile product offering**

- Multiple lengths, tip configurations, and diameters to meet a variety of clinical needs



## ORDERING INFORMATION

GUIDERIGHT™						
ORDER NUMBER	DESCRIPTION	DIAMETER (IN)	LENGTH (CM)	"J" RADIUS (MM)	FLEXIBLE TIP LENGTH (CM)	UNITS PER PACK
404555	Standard, Straight, "J"	0.035	150	3	5	5
404560	Super Stiff, "J"	0.032	180	3	8	5
404568	Standard, "J"	0.025	150	3	5	5
404569	Standard, "J"	0.032	150	3	5	5
404570	Standard, "J"	0.035	150	3	5	5
404571	Standard, "J"	0.035	180	3	5	5
404577	Standard, Exchange, "J"	0.035	260	3	5	5
404578	Firm, "J"	0.035	150	3	3	5
404610	Standard, "J"	0.038	150	3	5	5
404617	Standard, Exchange, "J"	0.038	260	3	5	5
404840	Standard, Straight	0.035	150	N/A	6	5
404843	Bentson, Straight	0.035	150	N/A	15	5
404844	Standard, Exchange, Straight	0.035	260	N/A	6	5
404878	Standard, Exchange, "J"	0.032	260	3	5	5

Caution: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available), at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. **Information contained herein for DISTRIBUTION in Europe, Middle East and Africa ONLY.** Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Tests performed by and data on file at Abbott. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

**Abbott Vascular International BVBA**

Park Lane, Culliganlaan 2B, 1831 Diegem, Belgium, Tel: 32.2.714.14.11

™ Indicates a trademark of the Abbott Group of Companies.

[www.cardiovascular.abbott](http://www.cardiovascular.abbott)

©2019 Abbott. All rights reserved. 2-EH-1-9263-01 03-2019

