

# INTRACARDIAC IMAGING

**INTRACARDIAC ECHOCARDIOGRAPHY (ICE)**

ViewMate™ Ultrasound Console .....E.1  
ViewFlex™ Xtra ICE Catheter ..... E.2

# ViewMate™

## Ultrasound Console

### Description

The ViewMate™ Ultrasound Console is a fully featured, compact console that quickly delivers high fidelity imaging to assist clinicians performing complex electrophysiology and interventional procedures. When the console is paired with a ViewFlex™ Xtra ICE Catheter, clinicians may visualize cardiac structures, blood flow and other devices within the heart.

### Product Highlights

- Fully connectable with DICOM<sup>®</sup> and network capabilities
- Multiple imaging modes, including 2-D, M-Mode, Pulsed Wave, Color Doppler and others
- Software-based architecture for rapid image optimization
- Enhanced far-field resolution enables assessment of distal devices and structures
- Up to three transducer ports allow for quick transducer selection and switching
- Easy-to-move console with small footprint

### Ordering Information

REORDER NUMBER	DESCRIPTION
H701336	ViewMate™ Ultrasound Console with Integrated Battery Pack
H702084	P8-3 TEE: Transesophageal Transducer
600023051	A2CW: Continuous Wave Transducer
H702488	L8-3: Linear Transducer
H701554	P4-1C: Transthoracic Transducer
600023347	ECG Accessory Cable
H701563	ViewFlex™ Catheter Interface Module (CIM) Mounting Kit, OUS
H701374	ViewFlex™ Catheter Interface Module (CIM) for ViewMate™



# ViewFlex™ Xtra

ICE Catheter

## Description

The ViewFlex™ Xtra ICE Catheter is a temporary intracardiac ultrasound catheter intended for use in adult and adolescent pediatric patients to visualize cardiac structures, blood flow and other devices within the heart when connected to a ViewMate™ ultrasound console via the ViewFlex™ Catheter Interface Module.

## Product Highlights

- Single-handed, self-locking steering
- 9 F catheter size
- 10 F introducer
- 90 cm insertable length
- Four-way tip deflection, anterior/posterior, left/right
- 120-degree tip deflection
- 64-element phased array

## Ordering Information

REORDER NUMBER	DESCRIPTION
D087031	ViewFlex™ Xtra ICE Catheter
H701374	ViewFlex™ PLUS ICE Catheter Interface Module for ViewMate™ Z Intracardiac and ViewMate™ Ultrasound Consoles
H701375	ViewFlex™ Catheter Interface Module for Philips‡ CX50 System

